



Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
FAPAS Food Chemistry Proficiency Test - 04520 - Multi - Mycotoxins in Maize Flour - July - September 2024	88.2	108	1.0
BIPEA 31B - Mycotoxins - Cereals - Triticale - 3 - 1931 - June 2024	33	27	-0.60
FAPAS Food Chemistry Proficiency Test - 04510 - Multi - Mycotoxins in Maize Flour - April - May 2024	97.1	120.94	1.1
BIPEA 31E - Mycotoxins - Feed - Animal Feed - 14 - 5 - January 2024	57	39	-1.13



BIO-SHIELD

FAPAS Food Chemistry Proficiency Test - 04500 - Multi - Mycotoxins in Durum Wheat - December 2023 - January 2024	131	127	-0.1
FAPAS Food Chemistry Proficiency Test - 22215 - Multi - Mycotoxins in Wheat Flour - November 2023 - January 2024	84.2	78.2	-0.3
BIPEA 31B - Mycotoxins - Cereals - Ground Corn - 6 - 187 - December 2023	496	357	-1.11
FAPAS Food Chemistry Proficiency Test - 04482 - Multi - Mycotoxins in Maize Flour - July - September 2023	186	158.7	-0.7
BIPEA 31C - 20 Mycotoxins - Baby Food - 21 - 3931 - June 2023	57	45	-0.75
BIPEA 31B - 185 - Mycotoxin Contamination in Cereals - Rice - 06-4331 - June 2023	62	50	-0.71
FAPAS Food Chemistry Proficiency Test - 04472 - Multi - Mycotoxins in Maize Flour April - June 2023	136.0	150.3	0.5
BIPEA 31B - 183 - Mycotoxin Contamination in Ground Barley - 05-0531 - February 2023	83.0	101.0	0.8
BIPEA 31E - 3 - Mycotoxin Contamination in Animal feed - 13-3031 - January 2023	238.0	241.0	0.05
FAPAS Food Chemistry Proficiency Test - 04461 - Mycotoxin Contamination in Animal Feed (Cereal Based) - December 2022 - January 2023	712.0	658.05	-0.5



FAPAS Food Chemistry Proficiency Test - 22199 - Multi - Mycotoxins in Wheat Flour - November 2022 - January 2023	81.9	71.86	-0.6
BIPEA 31B - 182 - Mycotoxin Contamination in Ground Corn - 18-0731 - December 2022	438.0	363	-0.68
FAPAS Food Chemistry Proficiency Test - 04448 - Multi - Mycotoxins in Maize Flour - July - September 2022	108.0	123.8	0.6
BIPEA 31B - 180 - Mycotoxin Contamination in Rice - 05 - 4331 - June 2022	50.0	39.0	-0.79
FAPAS Food Chemistry Proficiency Test - 04439 - Multi - Mycotoxins in Maize Flour - April - May 2022	94.3	101.2	0.3
FAPAS Food Chemistry Proficiency Test - 22187 - Multi - Mycotoxins in Wheat Flour - November 2021 - January 2022	74.4	74.21	0.0
FAPAS Food Chemistry Proficiency Test - 04426 - Multi - Mycotoxins in Maize Flour - October - December 2021	295.0	321.08	0.5
FAPAS Food Chemistry Proficiency Test - 04418 - Multi - Mycotoxins in Maize Flour - July - September 2021	122.0	132.75	0.4
FAPAS Food Chemistry Proficiency Test - 04411 - Multi - Mycotoxins in Maize Flour - April - June 2021	105.0	106.85	0.1



BIO-SHIELD

FAPAS Food Chemistry Proficiency Test - 22177 - Multi - mycotoxins in Wheat Flour - November 2020 - January 2021	75.6	79.31	0.2
FAPAS Food Chemistry Proficiency Test - 22175 - Deoxynivalenol, Zearalenone, T-2 and HT-2 in Animal Feed - September - October 2020	146.0	144.72	0.0
FAPAS Food Chemistry Proficiency Test - 04384 - Multi - Mycotoxins in Maize Flour - April - May 2020	87.3	96.3	0.5
FAPAS Food Chemistry Proficiency Test - 04378 - Multi - Mycotoxins in Maize Flour - February - March 2020	149.0	156.11	0.2
FAPAS Food Chemistry Proficiency Test - 04375 - Mycotoxin Contamination in Animal Feed (Cereal Based) - December 2019 - January 2020	906.0	854.0	-0.4
FAPAS Food Chemistry Proficiency Test - 22166 - Deoxynivalenol, Zearalenone, T-2 and HT-2 in Wheat Flour - November 2019 - January 2020	76.2	83.01	0.4