

# MILK & DAIRY ANTIBIOTICS



# A BINARY ANSWER IN THE SHORTEST POSSIBLE TIME



Test cold milk right from the truck, refrigerator or bulk tank

Option for visual interpretation

No need for extra disposables

## QUANTUM BT-CEF

W1030/W1060

Detection of  $\beta$ -lactams, including Cefalexin, and Tetracyclines in cow, sheep and goat milk

- Final Results in 5 min
- Validated by: ILVO | LILC




## QUANTUM BT-S

W1130/W1160

Detection of  $\beta$ -lactams, Tetracyclines and Sulfonamides in cow, sheep and goat milk

- Final Results in 5 min
- Validated by: ILVO



## QUANTUM BT-SQ

W1230/W1260

Detection of  $\beta$ -lactams, Tetracyclines, Sulfonamides and Quinolones in cow, sheep and goat milk

- Final Results in 5 min



## QUANTUM BT-SC

W1330/W1360

Detection of  $\beta$ -lactams, Tetracyclines, Streptomycin and Chloramphenicol in cow, sheep and goat milk

- Final Results in 5 min



## QUANTUM BL

W1430/W1460

Detection of  $\beta$ -lactams in fresh cow milk

- Final Results in 2 min



# FAST AND ACCURATE RAPID TESTS



Option for visual interpretation

Incubation at 40°C




## SYMMETRIC B+T

S1048/S1096

Detection of  $\beta$ -lactams and Tetracyclines in cow, goat and sheep milk

- Dipstick Format
- Validated by: ILVO | LILC | PIWET



## SYMMETRIC BTS

S1248/S1096

Detection of  $\beta$ -lactams, Tetracyclines and Sulfonamides in cow, goat and sheep milk

- Dipstick Format
- Validated by: ILVO





**FROM FARM TO CLOUD**


**Get results  
on the go  
with our  
Android  
Mobile App**



3 drops of cold milk



No incubation needed



Quantitative Results & GPS location

