

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Senf/Mustard R6152	Range for $z \leq 2$ according to R- Biopharm Ridascreen Fast Senf/Mustard R6152	Result (mg/kg)	Recovery (%)
FAPAS Celery & Mustard in Instant Soup Powder 27397A June - July 2024	43.4	21.7 - 65.1	44.0	101.4
FAPAS Celery & Mustard in Instant Soup Powder 27397B June - July 2024	NOT DETECTED	-	<1.0	-