

<b>Proficiency Test</b>	<b>Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Senf/Mustard R6152</b>	<b>Range for <math> z  \leq 2</math> according to R-Biopharm Ridascreen Fast Senf/Mustard R6152</b>	<b>Result (mg/kg)</b>	<b>Recovery (%)</b>
FAPAS Celery & Mustard in Instant Soup Powder 27397A June - July 2024	43.4	21.7 - 65.1	41.5	95.6
FAPAS Celery & Mustard in Instant Soup Powder 27397B June - July 2024	NOT DETECTED	-	<1.0	-
FAPAS Celery & Mustard in Instant Soup Powder 27361A June - July 2023	44.1	22.0 - 66.1	46.1	104.5
FAPAS Celery & Mustard in Instant Soup Powder 27361B June - July 2023	NOT DETECTED	-	<1.0	-
FAPAS Celery & Mustard in Instant Soup Powder 27328A June - August 2022	NOT DETECTED	-	<1.0	-
FAPAS Celery & Mustard in Instant Soup Powder 27328B	DETECTED	-	DETECTED	-



June - August 2022				
-----------------------	--	--	--	--

